

Versati™

Tabletop Centrifuge

Exceptional choice for reliable sample separation

Pelleting of large volume of temperature-sensitive heterogeneous mixtures are now easier, safer, and high-yielding with the Versati™ Tabletop Centrifuges. Esco tabletop centrifuge excel with its versatility, running features, and easy handling; developed with your molecular biology and IVF applications in mind.



Key Features



Maintenance-free, brushless motor



Superior safety features



Fast pre-cooling function



High temperature ramp rate
(for refrigerated models)



Large storage space (up to 99 programs)



Alarm system (10 different alarm tone)



Time-setting spin (up to 99 hours)

General Accessories for Tabletop Centrifuges



Swing-bucket Rotor

Aluminum swing-bucket rotor with circular flatbottom buckets made of polypropylene can hold up to 4 x 250 ml tubes. It has flexible adapters ideal for medical laboratories.



Microtiter Plate Rotor

This microtiter plate rotor has a maximum capacity of up to 6 plates. This can also accommodate deep well plate, culture plate, microtest/ terasaki plate, microsonic system and PCR well plate.



Aerosol-tight Fixed-angle Rotor

This TÜV Nord Certified Bioseal Rotor used for 1.5/2.0 ml tubes is also available in tabletop centrifuge models. Adapters are used to run 0.5 ml and 2.0 ml / 0.4 ml PCR tubes.



Fixed-angle Rotor

The maximum capacity of this fixed-angle rotor is 6 x 250 ml. It can also run tubes ranging from 1.5/2.0 ml to 50 ml using suitable adapters.

Note: There are a total of 12 rotor options for TCVI TCR and 47 available adapters for both models.

Also Available in: Micro Centrifuge Models



Micro Ventilated Centrifuge
(MCV-88-)



Micro Refrigerated Centrifuge
(MCR-88-)



Want to know more about this product? Scan the QR Code

Esco Micro Pte. Ltd.

21 Changi South Street 1 • Singapore 486 777 • Tel +65 6542 0833 • Fax +65 6542 6920
medical@escoglobal.com • esco-medical.com

Esco Global Offices: Bangladesh | China | Denmark | Germany | Hong Kong | Indonesia | Italy | Lithuania | Malaysia | Myanmar | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | UK | USA | Vietnam