

MIRI® II-12

Multiroom IVF Incubator

One Patient – One Chamber One Patient – One Chamber

The MIRI® II-12 multiroom IVF incubator dedicates one chamber for one patient to ultimately secure a minimum-stress environment for embryos. The chambers are completely independent of each other. Any disruption on one chamber (e.g. temperature drop after opening the lid) has zero impact on the rest of the system.



Key Features



Heated bottom

Provides direct heat transfer to the culture dishes for stable heat regulation.



Reliable gas mixing system

This allows gas phase flexibility. The gas mixer of MIRI® II-12 gives total control over CO₂ and O₂ concentration levels while also giving flexibility over what gas input is desired.



Fast Recovery

One important benefit is the speed of recovering temperature and gas parameters after opening a chamber.



Built-in pH System

The MIRI® II-12 has a built-in pH system that the user can use to calibrate the pH meter and measure the pH of the media.



Heated Lid

Prevents condensation and improves temperature regulation and recovery enhancing uniformity within each chamber.



Multiroom System

The MIRI® II-12 multiroom IVF incubator have multiple independent chambers with very stable environments, allowing embryologist to culture embryos from individual patients in individual chambers.



Stress-free validation of chamber parameters

The MIRI® II-12 can be connected to an external device such as the MIRI® GA for gas and temperature validation.



Touchscreen PC

Allows user to control and monitor important parameters (temperature, gas concentration) simultaneously and gives visual and audible alarm to signal critical conditions. It has large display that can be easily seen from a distance.



These limited-edition MIRI® II-12 units are available in line with Esco Medical's 10th year anniversary.

Want to know more about this product?

Scan the QR Code

www.esco-medical.com



Esco Medical, Aps

Kringelleg 10, 8250 Egaa, Denmark
medical@escolifesciences.com

Esco Micro, Pte. Ltd. (Headquarters)

19 Changi South Street 1, Singapore 486 779
Tel +65 6542 0833

Esco Global Offices: Bangladesh | China | Denmark | Germany | Hong Kong | Indonesia | Lithuania | Malaysia | Myanmar | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | UK | USA | Vietnam

9010451_MIRI 2 Flyer_A4_vF-011525

Esco can accept no responsibility for possible errors in catalogues, brochures and other printed materials. Esco reserves the right to alter its products and specifications without notice. All trademarks and logotypes in this material are the property of Esco and the respective companies.