

# MIRI<sup>®</sup> Laminar Flow Cabinet

A vertical laminar flow cabinet that provides a constant and reliable airflow giving superior quality on the air supplying the working zone. Thereupon provides sample protection thoroughly. The MIRI<sup>®</sup> Laminar Flow is equipped with a highly energy-efficient ECM blower which minimize power consumption without affecting the cabinet's performance.



## Key Features



### ISO Class 3

Equipped with HEPA/ULPA filter with efficiency of 99.999% towards micro particulate of 0.3 µm.



### Built-in Monitor

User may opt to have a touchscreen monitor equipped with a mini PC for better workflow.



### Energy Efficient

The latest generation ECM blower allows the cabinet to have minimal power consumption while maintaining the cabinet's performance.



### Microscope Provision

The microscope provision allows user to easily install a stereo-microscope unit at their convenience.



### Anti-Microbial Coating

The ISOCIDETM effectively inhibits and eliminates the growth of unwanted microbials at the surface of the cabinet's exterior.

## Optional Accessories



### Support Stand:

- With leveling feet
- With Caster Wheels
- Motorized stand with Caster Wheels



UV Kit

Want to know more about this product?

Scan the QR Code

[www.esco-medical.com](http://www.esco-medical.com)



### Esco Medical, Aps

Kringelled 10, 8250 Egaa, Denmark  
medical@escolifesciences.com

### Esco Micro, Pte. Ltd. (Headquarters)

19 Changi South Street 1, Singapore 486 779  
Tel +65 6542 0833

Esco Global Offices: Bangladesh | China | Denmark | Germany | Hong Kong | Indonesia | Lithuania | Malaysia | Myanmar | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | UK | USA | Vietnam

9010593\_Blue\_MLF Flyer\_A4\_vA\_121824

Esco can accept no responsibility for possible errors in catalogues, brochures and other printed materials. Esco reserves the right to alter its products and specifications without notice. All trademarks and logotypes in this material are the property of Esco and the respective companies.