

MIRI® Evidence

RFID Witnessing and Traceability System

The Ultimate Least Invasive RFID System.

The MIRI® Evidence RFID Witnessing and Traceability System streamlines scheduling and documentation in fertility clinics while helping prevent mismatches. It automatically arranges tasks at the start of treatment, updates schedules for any changes, and ensures accurate data registration. Missing data triggers automatic task creation for the responsible party.



Key Features



Hybrid Label

MIRI® Evidence ensures that the correct sample is always in use, as the machine-readable label is visually validated and cross-checked against the selected patient, effectively eliminating the risk of any mismatch.



Chain of Custody

MIRI® Evidence ensures precise tracking and documentation of every step, from procurement to recipient transfer. It automatically logs all equipment used and tracks events and movements, highlighting any protocol deviations in the “chain of custody” log.




Effortless Integration (Retrofit)

Designed with clinical flexibility in mind, the MIRI® Evidence easily integrates with existing equipment. Compatible with most workstations, our RFID pole reader ensures smooth operation across various laboratory setups.



Automatic Mismatch Prevention

The RFID technology system employed in the MIRI® Evidence will automatically detect any mismatch during a procedure once a labelled dish is placed within the area covered by the MIRI® Evidence Pole reader, Arc reader and Tube reader.



Least Invasive Solution for Your Workflow

MIRI® Evidence is designed to integrate seamlessly into your existing workflow. There's no need to alter your processes, adjust routines, or modify your tabletop setup. Everything remains the same—just more efficient and hassle-free.



Cycle Credit

The MIRI® Evidence concept ensures all containers for gametes and embryos are labeled, as labeling only key containers creates false safety.

Major: Ovum Pickups (IVF, ICSI, Transfer, and Cryo procedures)
Minor: Frozen Embryo Transfer + IU!

Want to know
more about
this product?

Scan the
QR Code

www.esco-medical.com



Esco Medical, Aps

Kringelld 10, 8250 Egaa, Denmark
medical@escolifesciences.com

Esco Micro, Pte. Ltd. (Headquarters)

19 Changi South Street 1, Singapore 486 779
Tel +65 6542 0833

Esco Global Offices: Bangladesh | China | Denmark | Germany | Hong Kong | Indonesia | Lithuania | Malaysia | Myanmar | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | UK | USA | Vietnam

9010595_Blue_MIRI Evidence_Flyer_A4_vD_051926

Esco can accept no responsibility for possible errors in catalogues, brochures and other printed materials. Esco reserves the right to alter its products and specifications without notice. All trademarks and logotypes in this material are the property of Esco and the respective companies.