

MIRI® M

Multiroom Incubator for IVF

This equipment is a CE-marked device and is in conformity with the essential requirements of the medical devices EU regulation 2017/745.

Key Features

Moveable Chambers

Maintains a stable $\pm 0.1^{\circ}\text{C}$ chamber temperature, even undocked for 30 mins.

Built-in Digital Traceability

Seamless logging of chamber movements for improved traceability and compliance.

Advanced Alarm System

MIRI® M incubator alerts users with pop-ups, icons, colors, and back panel backlights.

Practical User Interface

Integrated software enables easy monitoring and adjustment via a touchscreen menu.

Improved Multiroom System

- Gas recovery (~3 min for CO_2 , < 5 min for O_2)
- Temperature recovery (< 1 min)
- Temperature stability: $\pm 0.1^{\circ}\text{C}$
- Gas stability: $\pm 0.2\%$



Want to know
more about
this product?
Scan the QR Code!



Esco Medical, Aps
Kringelld 10, 8250 Egaa, Denmark
medical@escolifesciences.com

Esco Micro, Pte. Ltd. (Headquarters)
19 Changi South Street 1, Singapore 486 779
Tel +65 6542 0833



Improved Sample Transportation

Say goodbye to the risks of unstable handling and hello to a safer, more reliable way to transport embryos with confidence.

Engineered with an improved, state-of-the-art transportation system, MIRI® M Chambers provide users with a larger, more ergonomic surface to hold onto, significantly enhancing grip and control during manual carrying.